

Prontosan[®] BREAKS THE BIOFILM CYCLE

A proactive approach using a combination strategy of Prontosan Solution and Prontosan Gel X as part of wound bed preparation may prove helpful and aims to:

- Reduce the biofilm burden (Prontosan Solution)
- Prevent reconstitution of the biofilm (Prontosan Gel X)

Appropriate Time for Appropriate Wound DESCRIPTION OF WOUND OBJECTIVE HOW TO USE ACUTE WOUND - SURGICAL PRIMARY & SECONDARY INTENTION HEALING Rinse with Solution ■ High risk patient* Cleans ■ No slough Prevents biofilm/complications Minimal exudate ACUTE WOUND e.g. trauma Soak with Solution Debris Cleans Haematoma Prevents biofilm/complications CHRONIC WOUND - GRANULATING Soak with Solution Consider Gel X ■ High risk patient* Low exudate Prevents biofilm/complications





CHRONIC WOUND Soak with Solution Apply Gel X

- Light slough
- Low exudate



Cleans

Prevents biofilm / complications




CHRONIC WOUND - CRITICALLY COLONISED/INFECTED Soak with Solution

- Medium/high exudate
- Static wound
- Slough



- Prevents biofilm/complications





Apply Gel X

Compatible with other dressings | Once opened, Prontosan Solution and Gel X can be used for up to 8 weeks (single patient use)

^{*}High risk patient: Co-morbidities such as Diabetes, immuno-compromised, steroidal use, patients with previous wound infection and or biofilm and slough.