B. Braun Closure Technologies
Manipler®AZ – 35W

US	GB	Instructions for use
		Single use skin stapler
	DE	Gebrauchsanweisung
	FR	Einweg-Hautklammergerät
	FR	Mode d'emploi
		Agrafeuse cutanée à usage unique
	ES	Instrucciones de manejo
		Grapadora cutánea desechable
	IT	Istruzioni per l'uso
		Suturatrice cutanea monouso
	PT	Instruções de utilização
		Aparelho de clamp cutâneo descartável
	NL	Gebruiksaanwijzing
		Wegwerp-huidcliptapparaat
	NO	Bruksanvisning
		Hudstiftemaskin til engangsbruk
	SE	Bruksanvisning
		Engångs hudclipsapparat
	FI	Käyttöohjeet
		Kertakäyttöinen ihohakaslaitte
	HU	Használati utasítás
		Egyszer használatos sebkácsoló készülék
	PL	Instrukcja obsługi
		Stapler skórný jednorazowego użytku
	RU	Инструкция по применению
		Инструмент для наложения кожных скоб
	TR	Kullanım Talimatları
		Tek kullanımlık cilt stablerli
	RO	Instrucţiuni de folosire
		Dispozitiv de sutură de unică folosinţă
	DK	Bruugsvejledning
		Engangs hudstapler

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Maniplier®AZ–35W Single use skin stapler

Description

The Maniplier®AZ is a sterile, single use skin stapler designed to deliver 35 stainless steel staples for routine wound closure.

Indications

The Maniplier®AZ skin stapler can be used for routine skin closure in a wide variety of surgical procedures.

Applications

- Prior to using the skin stapler verify compatibility of skin stapler remover.
- Using sterile technique remove the skin stapler from the package.
- With a tissue forceps, pull skin edges together until edges evert (Fig.1)
- Position the skin stapler over the center of the everted skin edges (Fig.2)
- To form the staple, squeeze the trigger firmly together until the trigger motion is halted (Fig.3)
- Release the trigger (Fig.4)
- Move the stapler off the incision (Fig.5)
- Always verify the correct position of the staple.

Contra-Indications

- When it is not possible to maintain at least 5mm distance from the stapled skin to underlying bones, vessels, or internal organs, the use of staples for skin closure is contraindicated.
- The Maniplier®AZ staples contain nickel, which might lead to allergic reactions with some patients.
- Do not reuse or re-sterilize. This may cause damage on functional capacities of the product and product and may cause damage or infection in patients.

Adverse Reactions

Adverse effects associated with the use of this device include: wound dehiscence; allergic response in patients with known sensitivities to metals contained in 316L stainless steel (i.e., chromium, nickel and iron) enhanced bacterial infectivity; minimal acute inflammatory tissue reaction; and pain, edema, and erythema the wound site.

MRI Safety Information


The Staple of Maniplier®AZ is MR Conditional.

- A) Non-clinical testing demonstrated that the Staple of Maniplier®AZ is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:
- Static magnetic field of 3-Tesla or less
 - Maximum spatial gradient magnetic field of 1,500-Gauss/cm (extrapolated) or less
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence)
 - Normal Operating Mode of operation for the MR system

B) MRI-Related Heating

In non-clinical testing, the Staple of Maniplier®AZ produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in a 3-Tesla MR system (3- Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI):-

MR system reported, whole body averaged SAR	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.2-W/kg
Highest temperature change	1.8°C
Temperature scaled to whole body averaged SAR of 4-W/kg	2.5°C

C) Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Staple of Maniplier®AZ. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e. as seen on the gradient echo pulse sequence) extends approximately 15-mm relative to the size and shape of the Staple of Maniplier®AZ.

Pulse Sequence	TI-SE	TI-SE	GRE	GRE
Signal Void Size	271-mm ²	175-mm ²	651-mm ²	467-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Safety instructions

- Do not use the Maniplier®AZ if the date of expiry has passed.
- Store the Maniplier®AZ at room temperature and use it directly after the package has been opened.
- Never use products from open or damaged sterile packaging.
- Do not re-sterilize, as this might impair the functioning of the product and injure the patient.
- Use the Maniplier®AZ skin stapler only on one patient.
- After using, dispose of the product properly as medical waste.

Sterility

Sterilization is effective until the date of expiry, if the package is stored correctly in a closed and intact package.
Shelf life is 5 years.

Delivery formats and package sizes

The Maniplier®AZ –35W is delivered in a box of 6 units.

Explanation of the symbols on the packaging

	Single use
	Use by date
	Sterilized by ethylene oxide
	CE marking and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC
	Batch code
	See instructions for use
	Catalogue number Art no. 783100
	Staple size 6.9 mm x 4.2 mm

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Title: DMA/HOJ/MANI/MANIPLER AZ Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)
Title: SUT-PR-ES04-Técnico Sistemas Impresión
Date: Wednesday, 24 February 2021, 08:20 W. Europe Daylight Time
Meaning: Document signed as Author

UserName: García, Guillermo (garcgues)
Title: SU-PR-ES04-Technical support manager
Date: Wednesday, 24 February 2021, 08:25 W. Europe Daylight Time
Meaning: Precheck of Document

UserName: Alonso, Adriana (alonades)
Title:
Date: Thursday, 25 February 2021, 08:39 W. Europe Daylight Time
Meaning: Approve Document

UserName: Piera, Eva (piermaes)
Title: Senior Regulatory Affairs
Date: Thursday, 25 February 2021, 10:56 W. Europe Daylight Time
Meaning: Approve Document

UserName: Laporte, Martina (lapomaes)
Title: Director Quality & Regulatory Affairs CoE CT
Date: Thursday, 25 February 2021, 23:17 W. Europe Daylight Time
Meaning: Final Release of the Document