
Parameter	Condition
Static magnetic field strength (B0)	1.5 T, 3 T or 7 T
Maximum spatial field gradient	69 T/m (6 900 Gauss/cm)
RF excitation	Circularly Polarized (CP)

Non-clinical testing has demonstrated the aneurysm clips are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions. Failure to follow these conditions may result in injury.

YASARGIL® Aneurysm Clip

Parameter	Condition
RF transmit coil type	There are no transmit coil restrictions.
Maximum head SAR (normal operating mode)	■ 1.5 T: 1.0 W/kg ■ 3 T: 1.1 W/kg ■ 7 T: 2.0 W/kg
B1+ _{rms}	■ 1.5 T: 4.0 μT ■ 3 T: 2.1 μT
Scan duration	15 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR image artifact	The presence of this implant may produce an image artifact.

YASARGIL® Aneurysm Clip

- B1+_{rms} shall be used on all MR systems with this limitation parameter.
- Use SAR only on MR systems not providing B1+_{rms} limitation.
- The CranioFix2 titanium clamps are MR Conditional. Non-clinical testing demonstrated that the CranioFix2 titanium clamp is MR Conditional. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:
 - Static magnetic field of 3 Tesla or less
 - Maximum spatial gradient magnetic field of 720 Gauss/cm or less
 - Maximum MR system reported, whole body average specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning (i.e., per pulse sequence)
 - First Level Controlled Operating Mode for the MR system
 - Do not take the titanium clamp applicator into the MR environment. It is MR Unsafe.

CranioFix2
Cranial Fixation System

Pulse Sequence	Signal void size	Plane orientation
T1-SE	503 mm²	Parallel
T1-SE	596 mm²	Perpendicular
GRE	701 mm²	Parallel
GRE	563 mm²	Perpendicular

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the CranioFix2 titanium clamp, therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. In non-clinical testing, the image artifact caused by the device extends approximately 5 mm from the CranioFix2 titanium clamp when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

Artifact information
Under the scan conditions defined above, CranioFix2 titanium clamps are expected to produce a maximum temperature rise of +2.1 °C after 15 minutes of continuous scanning.

MRi-related heating
Under the scan conditions defined above, CranioFix2 titanium clamps are expected to produce a maximum temperature rise of +2.1 °C after 15 minutes of continuous scanning.

CranioFix2
Cranial Fixation System

Neuro-Patch® Dura Substitution Product

MRi examinations do not present an additional risk to implant wearers of Neuro-Patch®.


For more information visit our website: www.bbraun.com

For detailed patient information on the products listed in this implant information and for other languages of the implant information, please visit our website: www.bbraun.com/neurosurgicalpassport or use the QR code

AESCULAP® – a B. Braun brand

Aesculap AG
Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600
www.bbraun.com

Neuro Plating System Cranial Fixation System


Non-clinical testing has demonstrated the Neuro Plating System is MR Conditional. A patient with this device can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial magnetic field gradient of 3,000 gauss/cm (30T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1 W/kg and head SAR of 1 W/kg for landmarks above the shoulder
- Normal Operating Mode for gradient output

Under the scan conditions defined above, the Neuro Plating System is expected to produce a maximum temperature rise of less than 5.00 °C after 7 minutes of continuous scanning. A cool down period of 5 minutes is needed after each 7 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends approximately 2 mm from the Neuro Plating System when imaged with a gradient echo pulse sequence and 3.0 T MRI system.

Manufacturer acc. to MDD 93/42/EEC of the Neuro Plating System is:
Osteonic Co., Ltd. / 1004Ho, 38, Digital-ro 29-gil, Guro-gu, Seoul, Korea

Optional space for implant labels

AESCULAP®

Implant Information

Patient Name

Address

ZIP-Code, City

Date of Birth (DD/MM/YY)

Date of Implantation

Healthcare Institution

Patient Record

Signature of Surgeon / Stamp

The owner of this information has one or more of the following devices implanted (please tick):


- ☐ YASARGIL® Aneurysm Clip
- ☐ CranioFix2
- ☐ Neuro Plating System
- ☐ Neuro-Patch®

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