Prontosan® – Wound Gel X

Manufacturer
Fabricant

Medical device
Dispositif médical

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Contains a medicinal substance
Contient des substances dangereuses

en Prontosan® – Wound Gel X

Intended Purpose: For cleansing, moistening and decontamination of acute, chronic and infected dermal wounds and burns.

1. Introduction: Chronic skin wounds are often coated with slough, necrotic tissue and/or biofilm. These coatings are difficult to remove and lead to delayed wound healing. Therefore, proper wound cleansing is essential. The use of Prontosan® Wound Gel X provides long-lasting cleansing and decontamination of the wound bed between dressing changes. Acute wounds also require proper cleansing as they are generally contaminated with debris and microorganisms. These contaminants can interfere with the normal wound healing process and lead to complications such as infection. For acute wounds that require suturing, Prontosan® Wound Gel X should be applied after surgical intervention and suturing. Due to the unique combination of ingredients (i.e. the antimicrobial substance, polyhexanide, and the surfactant, betaine), Prontosan® Wound Gel X is ideal for the prevention of biofilm formation. Test results support the claim that Prontosan® Gel X is an effective barrier to reduce microbial penetration through the dressing. For optimal results, Prontosan® Wound Irrigation Solution (see separate product information) should be used for cleansing the wound and the skin area around the wound prior to treatment with Prontosan® Wound Gel X. Prontosan® Wound Gel X contains a hydrocolloid that supports the creation of an optimal moist environment for wound healing, thereby reducing the likelihood of scarring. The summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed) and can be found using basic-UDI-DI: 4039239000009232Q or the product name. <https://ec.europa.eu/tools/eudamed>. Or requested from the Manufacturer.

2. Product profile and areas for use: For cleansing, moistening and decontamination of:

- acute non-infected and infected wounds: traumatic wounds such as lacerations, crush injuries (if suturing is indicated Prontosan® Wound Gel X should be applied after surgical intervention).
- chronic non-infected and infected wounds (e.g. vascular ulcers, diabetic ulcers, pressure ulcers, including complex, recalcitrant wounds, cavity wounds, difficult-to-access wounds).
- postoperative wounds.
- thermal and non-thermal burns degree: I, IIa, IIb and III.
- fistulas and abscesses.

3. General use: To ensure optimal cleaning results with Prontosan® Wound Gel X, aseptic application is recommended whenever possible. Prontosan® Wound Gel X should be copiously applied to the wound bed. If necessary, gently disperse Prontosan® Wound Gel X with a sterile spatula. Cavities and pockets should be filled with Prontosan® Wound Gel X. Dressings, gauzes, compresses or other absorbent wound fillers can be moistened with Prontosan® Wound Gel X before the dressing is applied. Prontosan® Wound Gel X may remain on the wound until the next dressing change. Depending on the frequency of dressing changes, varying amounts of Prontosan® Wound Gel X are applied. The surface of the wound should be kept continuously moist to ensure adequate cleansing and decontamination. Coatings are gently released and removed with the next dressing change. Application should be conducted frequently enough for all coatings to be readily removed and for optimal wound bed preparation.

4. Intended users: The Prontosan® Wound Gel X is intended to be used by professional users or by instructed adults only. The intended patient population are adults and children.

5. Tissue tolerability and biocompatibility: Dermatologically tested and evaluated as non-irritating and well-tolerated; painless; no inhibition of granulation or epithelialisation.

6. Side effects: In very rare cases, there may be a mild burning sensation after application of Prontosan® Wound Gel X, but this usually dissipates after a few minutes. Prontosan® Wound Gel X can cause allergic reactions such as itching (urticaria) and rashes (exanthema). In very rare cases (less than 1 out of 10,000), anaphylactic shock has been reported with Prontosan® Wound Gel X products.

7. Contraindications: Prontosan® Wound Gel X should not be used:

- if the patient is known to be allergic or if it is suspected that the patient may be allergic to one of the ingredients of the product.
- on the central nervous system or the meninges.
- in the middle or inner ear.
- in the eyes. If Prontosan® Wound Gel X comes into contact with the eyes, flush the eyes with running water and seek medical advice.
- on hyaline cartilage and in aseptic joint surgery. If Prontosan® Wound Gel X does come into contact with aseptic cartilage, it should be immediately irrigated with Ringer's solution or normal saline.

8. Restrictions of use: Pregnancy and lactation period: There is no evidence of mutagenic or embryo toxicity associated with the ingredients of this product. Due to the lack of relevant clinical trials and clinical experience with pregnant and breast-feeding women, Prontosan® Wound Gel X should not be used.

9. General safety instructions: Do not use for infusion, injection or ingestion. Do not use in combination with cleansing soaps (anionic surfactants), ointments, oils, etc. These substances should be thoroughly removed from the wound before use as they can affect the cleansing performance of Prontosan® Wound Gel X. Do not use damaged tubes. Keep tubes away from direct sunlight. Keep out of reach of children. The use of Prontosan® Wound Gel X does not replace the need for systemic or other adequate infection treatment.

10. Summary / technical information: Prontosan® Wound Gel X is a preserved product that provides a sustained antimicrobial barrier and has a shelf life of 8 weeks after opening. The tube should be closed immediately after use to prevent contamination. The top of tube should be protected from contamination during use. Tubes that have come into direct contact with the wound or have become contaminated in another way should be discarded. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Composition: Purified Water, Glycerol, Hydroxyethylcellulose, Betaine Surfactant, 0.1 % Polyaminopropyl Biguanide (Polyhexanide)*.

Appearance and smell: clear, colourless and virtually odourless, aqueous gel. **Shelf life:** according to the expiry date, store at room temperature. **Originality:** sterile, originality seal. **Tubes:** for single patient use. **Manufacturing site:** Made in Switzerland.

*if used separately, can be considered to be a medicinal product.

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B. Braun Australia Pty Ltd
**B. Braun Medical AG**