# Aqua B. Braun

*Ecotainer® Neo:*
*500 ml; 1000 ml:*

*Ecolav® Neo:*
*500 ml:*

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## en Aqua B. Braun Instruction for use for patients or caregivers

*Instruction for use for professional users:*
*Please refer to B. Braun homepage:*
*https://ifu.bbraun.com/ Or (QR code)*

**1000 ml solution contain:**
Aqua B. Braun is stored in the bottle and can be used without additional devices.

This packaging unit of sale contains multiple containers of irrigation solution. Such containers must not be distributed individually unless an instruction for use is provided with each container.

**Ecotainer® Neo:**
- Irrigation bottle with screw-cap.
- Aqua B. Braun is stored in the bottle and can be used without additional devices.

**Ecolav® Neo:**
- Irrigation bottle with twist-off cap.
- Aqua B. Braun is stored in the bottle and can be used without additional devices.

**Aqua B. Braun:**
- Provides mechanical cleansing for sterile irrigation of body cavities, tissues, wounds, and burns, including during and after surgery.
- Helps maintain a clear surgical field by removing blood, debris, and tissue.
- Is used for rinsing medical instruments and equipment.
- Provides humidity to breathing gases through respiratory air humidifiers.
- Do not reuse the product.

**Sterility**
Sterile and endotoxin-free irrigation solution. Sterilization is performed by moist heat.

Only use Aqua B. Braun if the tamper evidence seal/ container is undamaged and the solution is clear.

**Intended purpose**
Aqua B. Braun is a ready for use, sterile, single use medical device, intended to be used for irrigation applications (irrigation solution).

**Indication**
Aqua B. Braun is used as irrigation solution or rinsing agent in patients with the need for:
- surgical interventions (intra- and postoperative irrigation)
- rinsing of urinary catheters
- rinsing of wounds and burns
- moistening wound tamponades, cloths and dressings
- stoma care
- mechanical rinsing of the eyes
- respiratory air humidifiers

**Patient population**
Aqua B. Braun can be used for all patients under consideration of the respective indications, contraindications and warnings for use.

**Intended User**
The products should only be used by trained professionals such as healthcare workers, licensed doctors, nurses, pharmacists and caregivers. After proper medical assessment, instruction and documentation by healthcare workers, patients and caregivers may also handle the products according to national guidelines.

**Contra-indications**
*Do not use Aqua B. Braun:*
- To inject or infuse
- When the ingress or absorption of larger amounts of the solution is expected

**Side effects**
- Using Aqua B. Braun can cause burning sensations, irritations, or discomfort at the application site.
- Displaced catheters or drainage tubes can lead to unintended irrigation or infiltration of body structures.
- Too much volume or pressure can damage tissues or cause leaks.
- Contamination due to improper handling may cause infection, which can lead to sepsis.
- Continuous or excessive irrigation can cause fluid overload, leading to symptoms like kidney problems, breakdown of blood cells, blood clotting problems, hypothermia, nausea, vomiting, headaches, feeling dizzy, confused, changes in your vision, electrolyte imbalances and plasma dilution, such as low natrium, difficulty breathing, changes in blood pressure such as an increase in blood pressure, lung or brain oedema in severe cases, or seizures. These symptoms are particularly likely to happen after transurethral resection of the prostate (TURP).
If signs of fluid overload occur, the intervention will be stopped and appropriate treatment will be given.

**Precautions**
- Do not drop the container. If dropped, discard it.
- If you are unsure about using the irrigation solution, consult your general practitioner or pharmacist.
- Aqua B. Braun should not be combined or used with medicinal products.
- In hospital settings only, Aqua B. Braun may be warmed up to not more than 40°C for maximum two weeks.

**Warning**
Before irrigation in surgery or post-operative setting, healthcare practitioners will treat problems with your body's fluid or acid levels.
During irrigation in surgery or post-operative setting, healthcare practitioners will monitor for excessive fluid absorption. If fluid overload is suspected, irrigation will be stopped, and the patient will be evaluated. Symptoms may appear later.
The following recommendations should be followed: time of surgery should be reduced to the minimum (for procedures with increased risk of fluid absorption into the circulation, pressure should be kept at a minimum, ensure adequate blood coagulation, sufficient monitoring of fluid balance.
Special precautions are to be taken when handling the irrigation solution.
Do not save opened containers; discard any leftover solution.
Do not use Aqua B. Braun after the expiry date on the label.
Re-using single-use devices can lead to contamination and functional issues, posing risks to patients and users.
Keep out of sight and reach of children.

**Operating instructions**
Ecotainer® Neo:
1. Grip container and open container by removing the screw-cap.
2. Pour solution in a sterile kidney bowl or address patient arm.

Ecolav® Neo:
1. Grip container and open container by removing the twist-off cap.
2. Pour solution in a sterile kidney bowl or address patient arm by applying slight pressure on the soft bottle.

**Duration of use**
- The amount and duration of use depend on the surgical intervention and subsequent recovery. The container does not need to be completely emptied.
- To avoid significant absorption during surgery, your healthcare professional must avoid prolonged or continuous use of the solution.
- In hospital post-operative settings, irrigation may last more than 60 minutes, but not more than 24 hours continuously.
- For respiratory air humidification / nebulization, Aqua B. Braun is intended for use up to 30 days

**Disposal**
Ask your pharmacist how to dispose of the solution and containers no longer required.

**Storage and handling conditions**
Please refer to the storage information on the packaging.
In hospitals, compatible disinfectants can be used on the surface of irrigation containers.

**Clinical benefits fr Performances**
The clinical benefit of the irrigation solution is related to the result of the associated therapy, when used as intended.

**Notice to the user**
All serious incidents related to this product must be reported to B. Braun and the competent authority in the country where the product is being used. Instructions for Use for professional users are available on the B. Braun homepage: https://ifu.bbraun.com/

**Date of issue and if appl. revision date**
Date of the last revision 2026-06-29.

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