
Aqua B. Braun

0075/12635097/0826

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


- 1000ml:

Ecolav® Neo:
- 500ml:

1000 ml solution contain
Water for Injection: 1000 ml

This packaging unit of sale contains multiple containers of irrigation solution. Such containers shall not be distributed individually unless an instruction for use is provided with each container.

Ecotainer® Neo:
- Irrigation bottle with screw-cap.
- Aqua B. Braun is stored in a bottle and can be applied without interfaces to other devices.

Ecolav® Neo:
- Irrigation bottle with twist-off cap.
- Aqua B. Braun is stored in a bottle and applied without interfaces to other devices.

Aqua B. Braun exerts a mechanical cleansing action for sterile irrigation of body cavities, tissues, wounds and burns including intra- and postoperative surgical interventions. The irrigation solution provides the means to maintain a clear surgical field through the removal of blood, debris and tissue during surgical interventions.

Furthermore, Aqua B. Braun exerts a mechanical action for eye irrigation purposes and for rinsing of (medical) instruments and equipment. Finally, Aqua B. Braun is utilized as a humidifying agent for the purpose of moistening wound tamponades, cloths and dressings and in respiratory air humidifiers.

Do not re-use.

Sterility
Sterile and endotoxin-free irrigation solution.
Ecotainer® Neo/Ecolav® Neo: Sterilization is performed by moist heat.

Do not use if package is damaged.
Visually inspect the sterile barrier system for breaches and integrity prior to use.

Intended purpose
Aqua B. Braun is a ready for use, sterile, single use medical device, intended to be used for irrigation applications (irrigation solution).

Indication
Aqua B. Braun is used as irrigation solution or moistening agent in patients in need of
- surgical interventions (intra- and postoperative irrigation)
- rinsing of urinary catheters
- irrigation of wounds and burns
- moistening wound tamponades, cloths and dressings
- stoma care
- mechanical irrigation of the eyes
- respiratory air humidifiers

Patient population
Aqua B. Braun can be used for all patients under consideration of the respective indications, contraindications, and warnings for use.

Intended User
The products should only be used by authorized persons e.g.: Healthcare professionals, licensed and/or certified physicians, nurses, pharmacists, caregivers (according to local regulations) who have been adequately trained within the educational background in this technique.

Upon medical assessment and adequate instruction, patients and/or caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and/or caregivers are instructed in the correct use of Aqua B. Braun and that the instruction of caregiver and/or patient is documented by the healthcare worker according to local documentation guidelines.

Contra-indications
- Do not use for intravascular infusion, injection and/or intrathecal irrigation.
- Do not use for irrigation when absorption of larger amounts of the solution is to be expected (e.g. during transurethral electroresection (TUR) procedures or hysteroscopy).

Side effects
Side effects are possible.

Generally, displaced catheters or drainage tubes can lead to irrigation or infiltration of unintended structures or cavities.

Moreover, excessive volume or pressure during irrigation of closed cavities may result in distension or disruption of tissues. This may also include reperforation and anastomotic leaks. Inadvertent contamination from careless technique may transmit infection, which may lead to sepsis. Irrigation procedures may further cause burning sensations, irritations or other discomfort at application site.

If used for continuous irrigation with excessive volume or pressure or allowed to dwell inside large body cavities or large wound surfaces, there is a risk of absorption of larger amounts of the solution, that could lead to states of hypervolaemia and their secondary manifestations and complications.

Absorption of larger volumes into the circulatory system may result in:
- Renal insufficiency (especially where there is existing renal damage)
- Haemolysis
- Clotting disorders caused by inhibition of the clotting factor activity or by lowering the coagulation factor concentration through haemodilution
- Hypothermia
- Nausea, vomiting
- Neurological problems: headache, sleepiness, agitation, confusion, which can develop into a potentially fatal coma or convulsions
- Visual disturbances
- Hyponatraemia, hypocalcaemia, sometimes pronounced, reduced osmotic pressure, reduced levels of haematocrit and serum protein
- Bradycardia
- Increased blood volume indicated by dyspnoea, variations in blood pressure or even acute pulmonary or cerebral oedema

All these symptoms are linked to an excessive inflow of irrigation liquid into the systemic circulation and indicate hyperhydration which can be fatal. They can appear during intervention (such as long-lasting interventions, particularly when a large venous sinus is opened) or in the decremental phase.

Often seen after transurethral resection of the prostate (TURP), the described symptoms and signs have been referred to as the TURP syndrome, although they have also been described after other urological or gynaecological surgical procedures.

If this diagnosis is suspected, these symptoms require the intervention to be stopped after careful haemostasis. The natraemia and the haematocrit should be rapidly measured and appropriate treatment for hyperhydration with haemodilution should begin – notably by restricting the hydric flow and administering a suitable diuretic.

Precautions
- Only use the solution if the tamper evidence seal/container is undamaged and the solution is clear and colourless. Do not drop the container. Otherwise discard container.
- Plastic containers can occasionally be damaged during transport from the manufacturer to the hospital, but also within the clinic. This can lead to contamination with microbial or fungiform growth of the solution. A careful visual inspection is therefore always required before using the solution. Particular attention must be paid to minimal damage of the sterility barrier (i.e. to the closure) or missing integrity of the tamper evidence seal due to possible risk of contamination.
- Aqua B. Braun is not intended to be combined or used with medicinal products. Healthcare professional facing use of Aqua B. Braun irrigation solution in conjunction with medicinal products must check the safety information of the additive as provided by the respective manufacturer.
- The container can be warmed up at 40°C for maximum two weeks and must not be warmed above 40°C. Temperature of the product must be checked before usage based on the clinical need. The equipment used for warming must allow to regulate the required temperature.

Warning
Throughout the irrigation procedure when there is a risk of absorption of larger amounts of the solution regular checks must be carried out for excessive inflow of irrigation liquid into the systemic circulation. This is particularly suggested by the appearance of any of the symptoms characteristic for fluid overload (for detail please refer to section side effects). If absorption of clinically relevant amounts of fluid is suspected, the fluid administration should be stopped, and the patient evaluated for possible adverse consequences.

The following recommendations should be followed:
- Time of surgery should be reduced to the minimum especially for those procedures at the increased risk of excessive inflow (normally no excess of 60 minutes)
- Pressure in the operating cavity should be kept at a minimum
- Ensure adequate haemostasis
- Sufficient monitoring of fluid balance

Symptoms of fluid overload might also become apparent with some delay after the intervention.

Disturbances of the acid-base balance and disturbances of the water electrolyte balance should be corrected before the use of Aqua B. Braun in irrigation procedures accessing body cavities or intraoperative or postoperative rinsing.

Strictly sterile precautions are to be taken when handling the irrigation solution.

Aqua B. Braun is intended for single use only and does not have to be completely emptied.

Opened containers must not be saved, quantities left over must be discarded.


Re-use of single-use devices creates a potential risk for of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

The solution must not be used after the expiry date given on the label.
Store out of reach of children.


B BRAUN

Operating Instructions

Ecotainer® Neo:
1. Grip container and open container with screw-cap.
2. Pour solution in a sterile container (e.g. kidney bowl) or address directly to the patient area.

Ecolav® Neo:
1. Grip container and open container by twisting off the cap without touching the outlet area.
2. Pour solution in a sterile container (e.g. kidney bowl) or address directly to patient area by applying slight pressure on the soft bottle.

Duration of use

The volumes used and the duration of use depend on the extent of the intervention and on the postoperative course. Thereby the container does not have to be completely emptied.

However, to prevent for clinically relevant absorption in operative settings, prolonged or continuous treatment with the solution must be avoided.

In post-operative irrigation settings, the irrigation solution is used in the hospital and may last more than 60 minutes, but it is normally not continuously used for more than 24 hours.

For respiratory air humidification / nebulization, Aqua is intended for use up to 30 days.

Disposal

Disposal according to local guidelines and / or clinical protocols.

Storage and handling conditions

Please refer to the storage information on the packaging.
Aqua B. Braun can be warmed up in a water bath immediately before use.

The following disinfectants are compatible with the outer surface of irrigation containers:
- 2-Propanol 70%
- 80% Ethanol
- Chlorhexidine 2% (e.g. Softa® Cloth)
- Ethanol 74% + Propan-2-ol 10% (e.g. Softasept® N)
- 1-Propanol 70%
- Sterile 70% denatured ethanol blended with deionised water
- Mixture of isopropyl alcohol with hydrogen peroxide H_2O_2 30%

Clinical benefits & Performances

The clinical benefit of the irrigation solution is related to the result of the associated therapy, when used as intended.

Notice to the user

All serious incidents related to this product must be reported to B. Braun and the competent authority in the country where the product is being used.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: https://eifu.bbraun.com/.

Date of issue and if appl. revision date

Date of the last revision 2026-06-29.

Meaning of Symbols




	Manufacturer
	Date of manufacture
	Use-by date
LOT	Batch code
REF	Catalogue number
eifu.bbraun.com	Consult instructions for use
	Do not re-use
STERILE	Sterile fluid path
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside
	Do not use if package is damaged
	Fragile, handle with care
	Keep away from sunlight
	Keep dry
	Temperature limit
	Non-pyrogenic
MD	Medical device
UDI	Unique device identifier
	Importer
STERILE	Sterilized using steam or dry heat
Not for Injection.	Do not use for injection.

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B. Braun Melsungen AG