




Original Perfusor™ Line

Coloured coded minimum volume
extension tubing for use with
Syringe drivers

Original Perfusor™ Line

Colour coded minimum
volume tubing.




Assign high potential
drugs by colour



Quickly identify
specific drugs
at a glance

Excellent
infusion start-up
characteristics
when used with
syringe drivers

Low priming volume



Does not contain DEHP or PVC

Our innovative minimum volume tubing features colour-coded lines, allowing you to assign specific drugs to designated colours. This simplifies line identification in complex therapies, enabling quick and clear recognition of what's being infused at a glance.

Features

- Luer Lock fitting enables compatibility to all syringe drivers
- Pressure resistant up to 2 bar
- Made with PE (Polyethylene), does not contain PVC, DEHP or LATEX
- Small tube diameter (1.0 x 2.0 mm) to reduce residual volume
- Available in 3 colours and two lengths
- Low sorbing properties
- Excellent start-up characteristics
- Kink resistant



POLYETHYLENE
(PE)



Not manufactured
with DEHP









Not manufactured
with PVC



2 BAR PRESSURE
RESISTANT

Ordering Information

Product	Colour	Material	Connection	Tube Length (cm)	Diameter (mm)	Priming volume (ml)	Units per box	Code no. (REF)
Perfusor™ Line	 Blue	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723020
	 Blue	PE (PVC-free)	Luer-Lock	200	1.0 x 2.0	1.67	100	8723021
	 Green	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723023
	 Green	PE (PVC-free)	Luer-Lock	200	1.0 x 2.0	1.67	100	8723024
	 Magenta	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723025
	 Magenta	PE (PVC-free)	Luer-Lock	200	1.0 x 2.0	1.67	100	8723026

B. Braun Australia Pty Ltd | Norwest Business Park NSW, 2153
Tel. 1800 251 705 | vetcare.au@bbraun.com | bbraun-vetcare.com.au