


Original Perfusor™ Line

Coloured coded minimum volume extension tubing for use with Syringe drivers

Features

- Luer Lock fitting enables compatibility to all syringe drivers
- Pressure resistant up to 2 bar
- Made with PE (Polyethylene), does not contain PVC, DEHP or LATEX
- Small tube diameter (1.0 x 2.0 mm) to reduce residual volume
- Available in 4 colours and two lengths
- Low sorbing properties
- Excellent start-up characteristics
- Kink resistant



Ordering Information

Product	Colour	Material	Connection	Tube Length (cm)	Diameter (mm)	Priming volume (ml)	Units per box	Code no. (REF)	UV protected
Perfusor™ Line	Blue	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723020	No
	Blue	PE (PVC-free)	Luer-Lock	200	1.0 x 2.0	1.67	100	8723021	No
	Green	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723023	No
	Green	PE (PVC-free)	Luer-Lock	200	1.0 x 2.0	1.67	100	8723024	No
	Magenta	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723025	No
	Magenta	PE (PVC-free)	Luer-Lock	200	1.0 x 2.0	1.67	100	8723026	No
	Orange	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723017	Yes