



ASKINA[®] CARBOSORB

CHARCOAL DRESSING

B. Braun Wound Care. Right. From the start.

Askina® Carbosorb

Charcoal dressing

Askina® Carbosorb is a conformable, sterile wound dressing with a thick, absorbent layer of non woven viscose-rayon, and an activated charcoal cloth layer for bacteria and malodorous absorption.

Askina® Carbosorb has limited absorption capacity: it can be used alone on low exuding malodorous wound, or in combination with other absorbing dressings for moderately to heavily exuding malodorous wounds.

- Odour uptake at 75% relative humidity: 0.60-0.75g/cm²
- Absorbency: 20-40g/100cm²

INDICATIONS

Askina® Carbosorb may be used for the management of low to heavily exuding, partial to full thickness chronic wounds, in particular when infection and/or offensive odour is present.

- Venous leg ulcers
- Arterial leg ulcers
- Pressure injuries
- Fungating wounds
- Infected and/or malodorous surgical or traumatic wounds



ADVANTAGES

- Can be cut to size
- Comfortable
- Can be used alone on low exuding malodorous wounds
- Can be combined with other absorbing dressings for moderate to heavily exuding malodorous wounds

Askina® Carbosorb dressings can be fixed with compression stockings for wounds on legs or with bandages or tapes.

Askina® Carbosorb can be easily cut for custom sizing to suit different shaped wounds for improved comfort.



B. Braun Ordering Information

Askina® Carbosorb

Size	Units/Box	Code
10 x 10 cm	10	9025006
10 x 20 cm	10	9025014