



ASKINA® SOFT | ASKINA® SOFT CLEAR

STERILE POST-OPERATIVE WOUND DRESSING

B. Braun Wound Care. Right. From the start.

Askina® Soft | Askina® Soft Clear

Sterile post-operative wound dressings



Askina® Soft is a sterile post-operative dressing combining a breathable fleece material and a non-adherent absorptive pad.
Absorption capacity: 1.89g


INDICATIONS


Management of post-operative wounds, and all sorts of acute wounds such as cuts, lacerations, sutures, abrasions.

- Elastic wound dressing
- Breathable fleece material
- Skin friendly adhesive
- Non-adherent wound contact layer
- Easy to apply, even on difficult locations

B. Braun Ordering Information

Askina® Soft

	Size	Pcs/Pack	Article number
	7.5 x 5 cm	50	9086480
	9 x 5 cm	50	9086501
	9 x 10 cm	50	9086510
	9 x 15 cm	40	9086528
	9 x 20 cm	30	9086536



Askina® Soft Clear is a sterile post-operative dressing combining a semi-permeable transparent film and a non-adherent absorptive pad.

- Moisture vapour transmission rate: 677g/m²/24hrs
- Spread of blood (to edge): 0.5 - 1.0g
- Spread of blood (to back): none


INDICATIONS

Management of post-operative wounds, and all sorts of acute wounds requiring a waterproof barrier and a visual control.

- Easy to apply, even on difficult locations
- Skin friendly acrylic adhesive
- Allows the patient to shower with the dressing in situ
- Bacteria proof barrier
- Allows visual control

B. Braun Ordering Information

Askina® Soft Clear

	Size	Pcs/Pack	Article number
	7.5 x 5 cm	50	9086803
	9 x 10 cm	30	9086821