




Aesculap[®] Retrieval Bag

Laparoscopic Surgery

Retrieval Bag

Art. No.	Description
	
EJ022SU	Retrieval bag with memory wire • Easy intra-abdominal opening and insertion of tissue • Simple intra-abdominal closure • Especially stable, transparent bag • Tapered bag design for easier removal • For use with 10 mm trocars • Insertion sleeve for easy insertion • Bag volume 210 ml

Art. No.	Description
	 Detachable retrieval bag • Self-opening • Bag can be opened and closed intraoperatively • Bag can be intraoperatively detached • Especially stable, transparent bag • For use with 10 mm trocars • Insertion sleeve for easy insertion
EJ023SU	 • Self-opening with metal ring • Bag volume 260 ml
EJ024SU	 • Self-opening with metal ring • Bag volume 720 ml

Art. No.	Description
----------	-------------



EJ025SU

"Simple" retrieval bag

- Simple, low-cost retrieval bag
- Red cord for secure bag holding
- Especially stable, transparent bag
- For use with 10 mm trocars
- Insertion sleeve for easy insertion
- Bag volume 200 ml





Art. No.	Description	Size	Volume	Trocar Ø	Pieces/Pack
EJ022SU	Retrieval bag with memory wire	85 x 160 mm	210 ml	10 mm	10
EJ023SU	Detachable retrieval bag	90 x 150 mm	260 ml	10 mm	5
EJ024SU	Detachable retrieval bag	130 x 170 mm	720 ml	10 mm	5
EJ025SU	"Simple" retrieval bag	76 x 150 mm	200 ml	10 mm	10

Sterile Single
use product

AESCULAP[®] – a B. Braun brand

B. Braun Australia Pty Ltd | Norwest Business Park NSW, 2153 | Tel. 1800 251 705 | info.au@bbraun.com | bbraun.com.au
B. Braun New Zealand | PO Box 37353, Parnell, Auckland 1151, New Zealand | Tel. 0800 227 286 | info.nz@bbraun.com | bbraun.co.nz

The main product trademark "Aesculap" is a registered trademark of Aesculap AG. Subject to technical changes.
All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products.
Reproduction in any form partial or otherwise is not permitted.