



Art. No.	Description	Size	Volume	Trocar Ø	Pieces/Pack
EJ022SU	Retrieval bag with memory wire	85 x 160 mm	210 ml	10 mm	10
EJ023SU	Detachable retrieval bag	90 x 150 mm	260 ml	10 mm	5
EJ024SU	Detachable retrieval bag	130 x 170 mm	720 ml	10 mm	5
EJ025SU	"Simple" retrieval bag	76 x 150 mm	200 ml	10 mm	10

Sterile Single use product








**Aesculap® Retrieval Bag**  
 Laparoscopic Surgery


## AESCULAP® – a B. Braun brand

B. Braun Australia Pty Ltd | Norwest Business Park NSW, 2153 | Tel. 1800 251 705 | [info.au@bbraun.com](mailto:info.au@bbraun.com) | [bbraun.com.au](http://bbraun.com.au)  
 B. Braun New Zealand Pty Ltd | Level 8, 139 Quay Street, Auckland, 1010, New Zealand | Tel. 0800 227 286 | [info.nz@bbraun.com](mailto:info.nz@bbraun.com) | [bbraun.co.nz](http://bbraun.co.nz)

# Retrieval Bag

Art. No.	Description
	
<p><b>EJ022SU</b></p>	<p><b>Retrieval bag with memory wire</b></p> <ul style="list-style-type: none"> <li>• Easy intra-abdominal opening and insertion of tissue</li> <li>• Simple intra-abdominal closure</li> <li>• Especially stable, transparent bag</li> <li>• Tapered bag design for easier removal</li> <li>• For use with 10 mm trocars</li> <li>• Insertion sleeve for easy insertion</li> <li>• Bag volume 210 ml</li> </ul>

Art. No.	Description
	<p><b>Detachable retrieval bag</b></p> <ul style="list-style-type: none"> <li>• Self-opening</li> <li>• Bag can be opened and closed intraoperatively</li> <li>• Bag can be intraoperatively detached</li> <li>• Especially stable, transparent bag</li> <li>• For use with 10 mm trocars</li> <li>• Insertion sleeve for easy insertion</li> </ul>
<p><b>EJ023SU</b></p>	 <ul style="list-style-type: none"> <li>• Self-opening with metal ring</li> <li>• Bag volume 260 ml</li> </ul>
<p><b>EJ024SU</b></p>	 <ul style="list-style-type: none"> <li>• Self-opening with metal ring</li> <li>• Bag volume 720 ml</li> </ul>

Art. No.	Description
<p><b>EJ025SU</b></p>	 <p><b>"Simple" retrieval bag</b></p> <ul style="list-style-type: none"> <li>• Simple, low-cost retrieval bag</li> <li>• Red cord for secure bag holding</li> <li>• Especially stable, transparent bag</li> <li>• For use with 10 mm trocars</li> <li>• Insertion sleeve for easy insertion</li> <li>• Bag volume 200 ml</li> </ul>
