


# Aesculap® Aicon® Basket System













The Aesculap® AICON® Basket System features a streamlined design and ergonomic handling.

- A wavelike basket bottom designed to help keep instruments upright during the packing process.
- A smooth, rounded structure with no feet, suitable for use with soft wrap and rigid container processes.
- An optional lid that allows visual inspection of contents while enabling handling via the basket's side handles.
- A stackable configuration, with or without a lid.
- Compatibility with the AICON® Container System as well as existing Aesculap® container systems.
- Interchangeable labelling tag options that integrate with the AICON® Container System, available in two sizes and ten colours.



## Tags

56 x 29 mm

with lettering <sup>3</sup> Item no.	without lettering Item no.	Colors
JJ540B	JJ540U	Red 
JJ541B	JJ541U	Blue 
JJ542B	JJ542U	Green 
JJ543B	JJ543U	Yellow 
JJ544B	JJ544U	Gray 
-	JJ545U	Black 
JJ546B	JJ546U	White 
JJ546B-BC <sup>5</sup>		White 
JJ546B-IC <sup>6</sup>		White 
JJ547B	JJ547U	Orange 
JJ548B	JJ548U	Purple 
JJ549B	JJ549U	Brown 
















Colours may differ from the illustration

- 3| **with lettering** – one, two or three lines, 5 mm font size, up to 13 characters per line
- 5| 1<sup>st</sup> line **barcode (type 128A)**, 2<sup>nd</sup> line identical content of barcode as plain writing, 3<sup>rd</sup> line arbitrary plain writing; all up to 13 characters per line
- 6| 1<sup>st</sup> line **Instacount** barcode up to 10 characters, 2<sup>nd</sup> line arbitrary plain writing up to 20 characters, 3<sup>rd</sup> line arbitrary plain writing up to 13 characters

## Tags

56 x 19 mm

with lettering <sup>4</sup> Item no.	without lettering Item no.	Colors
JJ560B	JJ560U	Red 
JJ561B	JJ561U	Blue 
JJ562B	JJ562U	Green 
JJ563B	JJ563U	Yellow 
JJ564B	JJ564U	Gray 
-	JJ565U	Black 
JJ566B	JJ566U	White 
JJ566B-BC <sup>7</sup>		White 
JJ566B-IC <sup>8</sup>		White 
JJ567B	JJ567U	Orange 
JJ568B	JJ568U	Purple 
JJ569B	JJ569U	Brown 



Colours may differ from the illustration

- 4| **with lettering** – one or two lines, 5 mm font size, up to 13 characters per line
- 7| 1<sup>st</sup> line **barcode (type 128A)**, 2<sup>nd</sup> line identical content of barcode as plain writing, all up to 13 characters per line
- 8| 1<sup>st</sup> line **Instacount** barcode up to 10 characters, 2<sup>nd</sup> line arbitrary plain writing up to 20 characters


## Size 1/1

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
Basket lid			
JJ167R	528 x 253 x 18		
Basket			
JJ161R		537 x 250 x 35	535 x 248 x 32
JJ162R		537 x 250 x 55	535 x 248 x 52
JJ163R		537 x 250 x 80	535 x 248 x 77
JJ164R		537 x 250 x 105	535 x 248 x 102




## Size 3/4

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
Basket lid			
JJ367R	394 x 253 x 18		
Basket			
JJ361R		403 x 250 x 35	401 x 248 x 32
JJ362R		403 x 250 x 55	401 x 248 x 52
JJ363R		403 x 250 x 80	401 x 248 x 77
JJ364R		403 x 250 x 105	401 x 248 x 102



## Size 1/2

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
Basket lid			
JJ267R	231 x 253 x 18		
Basket			
JJ261R		240 x 250 x 35	238 x 248 x 32
JJ262R		240 x 250 x 55	238 x 248 x 52
JJ263R		240 x 250 x 80	238 x 248 x 77
JJ264R		240 x 250 x 105	238 x 248 x 102



## Size DIN

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
Basket lid			
JJ467R	476 x 253 x 18		
Basket			
JJ461R		485 x 250 x 35	483 x 248 x 32
JJ462R		485 x 250 x 55	483 x 248 x 52
JJ463R		485 x 250 x 80	483 x 248 x 77
JJ464R		485 x 250 x 105	483 x 248 x 102

B. Braun Australia Pty Ltd | Level 2, 4 Burbank Place, Norwest NSW 2153, Australia | Tel.1800 251 705 | info.au@bbraun.com | www.bbraun.com.au

B. Braun New Zealand | Level 8, 139 Quay Street, Auckland, 1010, New Zealand | Customer Care 0800 227 286 | Fax (09) 373 5601 | www.bbraun.co.nz