


Dermaslide® Needle


A definitive choice for skin closure

DERMASLIDE® DGMP*: Features of the needle with a micro-engraved body

- Excellent tissue penetration¹
- High bending performance¹
- Sharper after multiple stitches compared to a classical cutting needle²
- Designed to improve the stability of the needle in the grasp of the needle-holder³



Dermaslide® – PENETRATION PERFORMANCE³



*DGMP NOMENCLATURE:
D: 3/8 circle needle
G: Golf ball structure
MP: Micropoint needle

REFERENCES
1. B. Braun Test Report 113_139/A Benchmark Cosmetic Needles: Ethicon-Covidien -S&R.
2. B. Braun R&D Testing Services Report (RDR-22-13110-FJB-CMY).
3. B. Braun R&D Testing Services Report (PTR/02/FF/1202/JLT-B).