

Syringes

Injekt® Two-Piece Standard Syringes



- Green plunger for easy recognition of filling volume
- Black graduations for ideal contrast and readability
- High barrel transparency simplifies inspection of the content
- Polypropylene syringe barrel and plunger
- Latex and DEHP Free

Description	Size	Box Qty	Code
Luer Slip	1 mL	100	9166017V
Luer Slip	2 mL	100	4606027V
Luer Slip	5 mL	100	4606051V
Luer Slip	10 mL	100	4606108V
Luer Slip	20 mL	100	4606205V
Luer Lock	2 mL	100	4606701V
Luer Lock	5 mL	100	4606710V
Luer Lock	10 mL	100	4606728V
Luer Lock	20 mL	100	4606736V

Omnifix® Three-Piece Standard Syringes



- Plunger stopper with double sealing ring
- Easily noticeable backstop to prevent unintentional withdrawal of plunger from barrel
- Colour-coded dispenser boxes for quick identification by size
- Low stiction properties
- Polypropylene barrel and plunger
- Latex and DEHP free

Description	Size	Box Qty	Code
Luer Slip	1 mL	100	9161406V
Luer Slip	3 mL	100	4616025V
Luer Slip	5 mL	100	4616057V
Luer Slip	10 mL	100	4616103V
Luer Slip	20 mL	100	4616200V
Luer Slip	30 mL	100	4616308F
Luer Slip	50 mL	100	4616502F
Luer Slip	50/60 mL	100	4613503F
Luer Lock	3 mL	100	4617022V
Luer Lock	5 mL	100	4617053V
Luer Lock	10 mL	100	4617100V
Luer Lock	20 mL	100	4617207V
Luer Lock	30 mL	100	4617304F
Luer Lock	50/60 mL	100	4617509F

Omnican® N Insulin Syringes



- Highly transparent barrel with high contrast black graduation for easy reading and exact dosing
- Highly transparent barrel to make sure that air bubbles are easily identified
- Integrated thin-wall needle to support smooth injection

Description	Size	Box Qty	Code
Omnican® N 50	30G x 8 mm	100	FT9172210S
Omnican® N 100	30G x 8 mm		FT9172220S
Omnican® N 100	30G x 12 mm		FT9172230S
Omnican® N 40	30G x 12 mm		FT9172240S